



COMBI LOCK STOPPER CLOSE OFF DEVICE FOR LUER CONNECTORS



The Combi Stopper is suitable for secure closure of male and female Luer and Luer-lock connectors.

With the Combi Stopper, pre-filled syringes, infusion ports, 3 Way Taps, Extension Sets, IV sets, Needle Free Injection Sites, Manifolds etc for example, can be closed securely. They can also be used with Catheter Tip Syringes.

Available in three Colours, White Blue & Red they are ISO compatible with all Luer Fittings including Syringes.



Combi Stopper made from polyethylene
Meets ISO for Luer Connections
For closure of male and female Luer connectors
Colour: WHITE, BLUE OR RED
Individually sterile packed
TGA Listed ARTG 223994

SOLMED PTY LIMITED. Unit 10, 42-44 Abel St, Jamisontown, Penrith NSW 2750

PHONE: 02-47210371