



WHEELCHAIRS

Features:

- Light Weight Aluminium - 17.2kg
- 130kg Maximum user weight
- Quick release wheels
- Solid Tyres
- 18" seat width, 17" seat depth
- Swing in/out leg rests
- Swing back removable armrests
- Height adjustable back upholstery
- Adjustable castors
- Multi position rear wheels

Specifications:

- Dimensions:
- Seat Width: 45cm
- Seat Depth: 43cm
- Overall Width: 62cm
- Overall Height: 98cm
- Safe Working Load: 130kg
- Weight of Chair: 17.5kg
- Weight without wheels: 14kgLink

Link Wheelchair



SOLMED PTY LIMITED ABN 16 132 179 250

**Unit 10, 42-44 Abel Street, Jamisontown,
Penrith. NSW 2750 PH: 02-47210371**

sales@solmed.com.au www.solmed.com.au