



INTRAVENOUS INFUSION ADMINISTRATION SETS



THESE SETS ARE SUITABLE FOR FLUIDS & DRUGS (NOT BLOOD)

Latex Free Needle Access Site, 220cm, 20 Drops/mL Drip Chamber 15 Micron Filter, Vented Perforator with PALL Hydrophilic Filter, Roller Regulator and Kink Resistant Tubing.



Latex Free Needle Free Swabable Access Site, 220cm, 20 Drops/mL Drip Chamber, Vented Perforator with PALL Hydrophilic Filter, Roller Regulator and Kink Resistant Tubing.



Latex Free In-Line Burette with Vented Perforator, Needle Free Swabable Access Site, Roller Regulator and Click Clamp, Safety Float Shut Off Valve, 150mL Burette Chamber 25cm Tubing.



Latex Free Administration Burette with Vented Perforator, Needle Free Swabable Access Sites, Roller Regulator and Click Clamp, Safety Float Shut Off Valve, 150mL Burette Chamber 255cm Tubing, 60 Drops/mL Micro Drip Chamber 15 Micron.



SOLMED PTY LIMITED. Unit 10, 42-44 Abel St, Jamisontown, Penrith NSW 2750

PHONE: 02-47210371