



WHEELCHAIRS

Features:

- Ultra-light 11kg frame
- Folding backrest
- Attendant hand brakes
- Rear wheel locks
- 19" seat width
- Safety lap belt
- Swing-away, removable leg rests
- Height adjustable, flip-up footplates
- Padded armrests
- 100kg safe working load

Specifications:

Weight of Chair: 11kg frame (12.5kg with leg rests);
Overall Seat Width: 482.5mm;
Seat Width (Between Arms): 482.5mm;
Seat Depth: 407mm;
Backrest Height: 457mm;
Leg rest Height Range: 356-483mm (add cushion for additional height);
Overall Width (Footprint): 559mm;
Overall Depth (Footprint): 635mm;
Armrest Height (Seat to Top of Handle): 254mm;
Wheels: Front: 203 x 25.4mm Solid Castors;
Rear: 305 x 57mm;
Safe Working Load: 100kg;
One year limited warranty

Shopper 12 Transit



SOLMED PTY LIMITED ABN 16 132 179 250

Unit 10, 42-44 Abel Street, Jamisontown,
Penrith. NSW 2750 PH: 02-47210371

sales@solmed.com.au www.solmed.com.au