



INTRAVENOUS TRANSFUSION ADMINISTRATION SETS



THESE SETS ARE SUITABLE FOR FLUIDS & DRUGS & BLOOD)

Latex Free Needle Free Swabable Access Site, 220cm, 20 Drops/mL Drip Chamber 200micron, Vented Perforator with PALL Hydrophilic Filter, Roller Regulator and Kink Resistant Tubing.



Latex Free, 2 X Needle Free Swabable Access Sites, Non Return Valve (Anti-Reflux) 220cm, 20 Drops/mL Drip Chamber 200Micron, Vented Perforator with PALL Hydrophilic Filter, Roller Regulator and Kink Resistant Tubing and Click Shut Off Valve.



**Rapid Flow Manual Pump Transfusion Set
Latex Free, Needle Free Swabable Access Site, Roller Regulator and Click Clamp, 220cm, Flexible Soft Drip Chamber 200Micron, 20 Drops/mL**



SOLMED PTY LIMITED. Unit 10, 42-44 Abel St, Jamisontown, Penrith NSW 2750

PHONE: 02-47210371